# Instructions for use

EN	Instructions for use	2
DE	Einbauanleitung	5
FR	Guide de fabrication	8
ES	Instrucciones de uso	11
IT	Istruzioni per l'uso	14
DA	Brugsanvisning	17
SV	Bruksanvisning	20
NL	Gebruiksaanwijzing voor de vervaardiging	23
PT	Instruções de Fabrico	26

**Icex<sup>®</sup> Ergo**

*Life Without Limitations*

**ENGLISH** – Caution: OSSUR products and components are designed and tested to ISO 10328. Compatibility and compliance with this standard is achieved only when OSSUR products and components are used with other recommended OSSUR or authorized components.

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ENGLISH

## Dear Customer,

Thank you for choosing Icex Ergo as the socket system for your client. Please read the following carefully before you start.

Icex Ergo is designed along with the Icecast Anatomy casting device; together they are used to create a total surface bearing socket directly onto the residual limb of transtibial amputees. This streamlined procedure bypasses several stages common in traditional socket production, allowing for a considerable reduction in the delivery time of the prosthesis.

The Icex Ergo braid is packed in a sealed foil pouch, care should be taken to avoid puncturing the pouch as this will cause the braid to set prematurely. The Icex Ergo kit should be stored flat in a cool, dry place.

The Icex Ergo braid has a shelf life of 12 months from date of manufacturing; the „best before“ date is marked on the foil pouch. Before immersing in water, always make sure that the braid feels soft and pliant. Protective gloves should be worn at all times when handling the Icex Ergo braid as the resin will adhere firmly to unprotected skin and clothing; Cleaning the skin with alcohol immediately may help remove uncured resin.

The Icex Ergo socket attachment is designed to house the Icelock 600 series and can be used with ratchet, lanyard and valve applications. The socket attachment is available with both a male pyramid and 4-

hole attachments, both tested to 100kg (220Lbs).

Icex Ergo is intended for certified prosthetists that have attended a specific training seminar. The 18 months warranty and responsibility of Ossur is limited for the intended use of the Icex Ergo product and applies only if the recommended manufacturing procedures are followed precisely.

## Instructions

- Examine the Residual Limb and collect all relevant measurements.
- Select appropriate Icross Liner type and size. Refer to Icross selection guidelines and sizing instructions.
- Have user don Icross Liner ensuring that it is centred distally and has no trapped air.
- Prior to socket fabrication it is recommended that a pre-cast pressure assessment be carried out.

## Pre Cast Pressure Assessment

- Select and apply the appropriate Silicone Disc from the Icecast Anatomy casting kit.
- Screw the Icex Ergo Casting Pin onto the Icross.
- Screw the Icex Ergo Stabilising Tool onto the Casting Pin.

- Prepare Icecast Anatomy with Contour Cell™ shaping chambers deflated and Main Chamber at 20mmHg.
- Connect and roll the Icecast Anatomy over the residual limb.
- During rolling the pressure in the Main Chamber will increase, once in place reduce the Main Chamber pressure back down to 20mmHg.
- Inflate the Contour Cell shaping chambers the desired pressure (See Icecast Anatomy Pressure Chart for recommended pressures).
- Inflate the Main Chamber to the desired pressure (See Icecast Anatomy Pressure Chart for recommended pressures).
- Leave in place for 3 minutes as this will help reduce the residual limb to its minimal volume whilst also indicating if the casting pressures are comfortable for the amputee.
- Deflate the Contour Cell shaping chambers and remove the Icecast Anatomy from the residual limb and unscrew the Stabilizing Tool.
- Place Steel Ring in the Ring Assembly Tool.
- Don latex gloves.
- Open Icx Ergo foil pouch. Remove all packaging exposing the carbon braid.
- Immerse the braid in the water, gently squeeze to assist water penetration.
- The Icx Ergo braid will remain workable for up to 120 seconds.
- Press the distal socket connector firmly down into the Steel Ring ensuring that the clips engage over the distal socket connector.
- Position Icx Ergo braid over the pin and rotate to ensure appropriate alignment of Distal Socket Connector, for convenience the distal socket connector is fully visible during the socket assembly.
- Hold the collar of the Stabilizing Tool while firmly tightening down onto the Distal Socket Connector of the Icx.
- Roll Icx Ergo braids upwards over the residual limb and trim around the patella to enable easy removal.
- Remove wrinkles by pulling and moulding the braids in a proximal direction. The knee should be extended but relaxed.

### **Icx Ergo Fabrication**

- Apply the Icx Ergo Inner Cover, reverse it so the plastic coated side is exposed, position the seam to the posterior and pull it over the residual limb. Ensure it is wrinkle free before trimming to length 5mm beyond the Silicone Seal distally.
- Select appropriate Pressure Pads from Icx Ergo Pressure Pad set remove backing and adhere to Icx Inner Cover over the bony prominences. (Alternative Pressure Pads are available if required).
- Prepare Icecast Anatomy with the Contour Cell shaping chambers deflated and Main Chamber at 20mmHg.
- Prepare container of water at 30°C (85°F). Ensure it is deep enough to completely submerge the braid.
- Connect and roll the Icecast Anatomy bladder over the Icx Ergo braid ensuring the anterior guide is positioned up the anterior of the residual limb.
- During rolling the pressure in the Main Chamber will increase, once in place reduce the Main Chamber pressure back down to 20mmHg.
- Inflate the Contour Cell shaping chambers to the pre-selected pressures.
- Finally inflate the Main Chamber to the pre-selected pressure.
- Disconnect the pump handle and support the Casting Bladder on a stool.

- The Braids should be left to set for 10 minutes.
- Remove Icecast Anatomy and unscrew the Stabilising Tool.
- Mark the anterior and the distal tip of the patella on the Icx Ergo Braid.
- The Icx Ergo may now be removed and trimmed (See technical manual for further detail).
- Assemble the Icelock 600 locking mechanism according to the product instructions supplied.
- Full weight bearing and dynamic alignment may be performed 45 minutes after initial immersion in water.
- Once the socket has been fit to the amputee and the final trimlines established smooth the edges with fine sandpaper and apply the Thermoformable Edge Tape.

### Packing List

Product Name	Product Number	Pieces
Icx Ergo Braid with Pyramid or	X-961550	
Icx Ergo Braid with 4-hole connector	X-962550	1
Icx Ergo Pressure Pad Set	X-190150	1
Icx Inner Cover	K-190300	1
Icx Steel Ring	X-990610	1
Icx Thermoformable Edge Tape	K-190452	1
Icx Interface Sheath	K-190200	1
Icx Ergo Instructions	IFU 0180	
Latex Gloves		1 Pair

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