







# UNLOADER ONE. THE ORIGINAL. THE STANDARD. THE UNLOADER.

# 3 cheers for 3 points.

The efficacy of our 3-Point Leverage System has been clinically proven by multiple peer-reviewed studies. Unloader One takes our patented GII® technology to the next level.



# Rachet it up... or down.

Featuring a tool-free ratchet, Unloader One enables on-the-fly adjustable and repeatable unloading, helping patients better manage their unicompartmental OA knee pain.

# Barely there.

As the lightest and lowest-profile OA knee brace on the market, Unloader One is barely noticeable to the user and fits inconspicuously and comfortably under clothing.





# Migration is for the birds.

Other OA knee braces often migrate. Össur's patented Sensil® silicone liners are specially designed to comfortably secure the Unloader One in place, maximizing effectiveness.

30 S GUARANTEE

Prescribe the Unloader One for any patient. If they are not happy with it for whatever reason, they can return it within 30 days for a full refund.

# **CLINICALLY PROVEN OSTEOARTHRITIS PAIN RELIEF**



RECENT CLINICAL RESEARCH BY THE STEADMAN PHILIPPON RESEARCH INSTITUTE\* ESTABLISHED THAT THE UNLOADER ONE® DECREASES PAIN, IMPROVES FUNCTION AND REDUCES THE USAGE OF PAIN MEDICATION.

# STEADMAN PHILIPPON RESEARCH RESULTS

Patients reported a significant improvement in disability (per WOMAC score):


Knee compartment unloaded	WOMAC Score improvement			
Medial	16 points			
Lateral	8 points			

At 6 months, a significant percentage of patients reported a reduction in medication use:

% of patients reporting medication reduction	Type of medication
27%	Pain medication
31%	OTC anti-inflammatories
35%	Prescription anti-inflammatories

Median patient satisfaction: 8 (out of 10)

# **EASY TO FIT AND WEAR**



Alligator clips make strapping reversible for external or internal fastening



Pre-attached popliteal pad improves patient comfort and maintains proper strap positioning



Color-coded buckles for easy and intuitive brace application



Matching dots on straps enable easier re-assembly of brace after washing



<sup>\*</sup> Research available at www.ossur.com/unloader

# **AVAILABLE IN PLUS SIZES, JUST LIKE YOUR PATIENTS**



## Comfortable containment.

For a variety of reasons, many sufferers of OA knee pain are overweight. For these patients, containing excess soft tissue is essential if they are to be braced effectively and comfortably. That's why our new Unloader One Plus model is equipped with a wider thigh strap¹ than the standard version, plus a gastroc strap that helps secure the brace in place?

## Extra large support.

Bigger legs require bigger support. So, for the new Unloader One Plus model, we included a sturdier upright with an Adjustable Dynamic Joint that provides clinicians with additional leverage for optimal pain management. Multiple, peer-reviewed clinical and biomechanical studies have proven the efficacy of this technology.



TS Custon

# **UNLOADER ONE® ORDERING INFORMATION**

## **UNLOADER ONE® CUSTOM**

Color	Right	Left		
Grey	B-115501711	B-115601711		
Black + 12 Colors (see below)	B-115500111	B-115600111		

## **UNLOADER ONE® OTS**

Size	Leg Circumference*	Medial	Lateral		
S	12.2-14" (31-35.6cm)	B-240 <b>X</b> 19712	B-240 <b>X</b> 29712		
M	14-16.75" (35.6-42.5cm)	B-240 <b>X</b> 19713	B-240 <b>X</b> 29713		
L	16.75-20" (42.5-50.8cm)	B-240 <b>X</b> 19714	B-240 <b>X</b> 29714		
XL	20-24" (50.8-61cm)	B-240 <b>X</b> 19715 B-240 <b>X</b> 29715			
	* 6" below mid-patella	X = 5 (Right), 6 (Left)			

## **Short Version**

Size	Leg Circumference*	Medial	Lateral		
XS	9.8-12.2" (24.9-31cm)	B-241 <b>X</b> 19711	B-241 <b>X</b> 29711		
S	12.2-14" (31-35.6cm)	B-241 <b>X</b> 19712	B-241 <b>X</b> 29712		
M	14-16.75" (35.6-42.5cm)	B-241 <b>X</b> 19713	B-241 <b>X</b> 29713		
	* 6" below mid-patella	X = 5 (Right), 6 (Left)			

## Notes:

Accessories available: Straps, pads, liners, AMS wrap, buckles, alligator tabs Short Length: 14.5" (36.8cm); Standard Length: 15.5" (39.4cm); Long Length: 16.5" (41.9cm) No weight limit. Sizing based on thigh/calf circumference.

## **UNLOADER ONE® PLUS CUSTOM**

Color	Right	Left		
Grey	B-119501700	B-119601700		
Black + 12 Colors (see below)	B-119500100	B-119600100		

## **UNLOADER ONE® PLUS OTS**

Size	Thigh Circ*	Calf Circ*	Part Number		
L	21-25" (53-64cm)	14-16.75" (36-43cm)	B-245 <b>XY</b> 0004		
XL	25-29" (64-74cm)	16.75-20" (43-51cm)	B-245 <b>XY</b> 0005		
	* 6" above/b	X = 5 (Right) or 6 (Left) Y = 1 (Medial) or 2 (Lateral)			

## **Long Version**

Size	Thigh Circ*	Part Number	
L	21-25" (53-64cm)	14-16.75" (36-43cm)	B-244 <b>XY</b> 0004
XL	25-29" (64-74cm)	16.75-20" (43-51cm)	B-244 <b>XY</b> 0005
	* 6" above/b	X = 5 (Right) or 6 (Left) Y = 1 (Medial) or 2 (Lateral)	

## **Custom Colors**

SB	OB NB	GR	СН	RD	PK	WH	SV	сс	YW	вк
----	-------	----	----	----	----	----	----	----	----	----

 $SB = Sky \ Blue, \ OB = Ocean \ Blue, \ NB = Navy \ Blue, \ GR = Green, \ CH = Champagne, \\ RD = Red, \ PK = Pink, \ WH = White, \ SV = Silver, \ CC = Charcoal, \ YW = Yellow, \ BK = Black$ 

USA 800 233 6263
CANADA 800 663 5982
WEB WWW.OSSUR.COM

FOLLOW ÖSSUR ON











© 2011 Össur P-530226, 7.11